

NEW

**Deodorization and moisture absorbing mat
for X-ray protection apron**

MediCore

ProteClean®

with hanging it on!



Sweat will cause of bad
odor and mold of
X-ray protection apron

Specifications

Size: 60 x 37 x 4cm (when hung)
Mat: 0.7cm thick
Cover fabric: 100% polyester
Contains: 100% silica gel
Made in Japan


Moisture absorbing and deodorization by silica gel of SilicaClean®

We adopted SilicaClean® which sandwiches high-performance silica gel.



● Moisture absorbing sensor

The color of the sensor changes according to humidity. If the color of the sensor changes from Blue to Pink, it means that moisture has accumulated, and please dry the ProteClean.



● Shoulder Mat

The moveable mat takes care of your shoulder you feel uncomfortable about sweat.

● Front and back Mat

Two-sided structure of SilicaClean® fits the inside of X-ray protection apron and enhances moisture absorbing and deodorization effect.

● Spacer

ProteClean keeps air permeability with inside spacer.

MediCore